Loughborough University  The Centre for Biological  Engineering	Document Change Note					
Document Ref: HTA-QS-FORM/003	Version N°:	1.0	Issue Date:			

DCN No: CBE-HTA/DCN/00013

Document Reference(s):	Doc	Current Version No:1.0					
HTA-MI- SOP0010	Risk Management Contingency Planning						
Reason for and overview of the change requested:	Addition of information about temperature Monitoring system 'Koolzone'.						
	Update of SOP to include details of new Koolzone Temperature Monitoring System which continually monitors Cryostorage, Fridge & freezer temperatures. This data can be viewed in specific software which is monitored by designated personnel. The system sends E-mail & Text alerts to designated personnel when the temperature is out of set ranges so they can act accordingly.						
Revision Category:	All personnel performing this procedure/process must read this version & complete the necessary training						
	$\boxtimes$	All personnel performing this procedure/process must read this version					
	No significant changes to document content – no requirement for personnel performing this procedure/process to read this version or complete training						
Systems effected by change:	Nor	ne	Other Changes required to implement:	None			
Date: 6 <sup>th</sup> December 2021	Rais	eed by: C.Kavanagh	Date: 13/12/2021	Approved by : R. Thomas ( Person Designa  K. Coopman ( Designated In			